Restylane[®]

REFYNE[™]

Injectable Gel with Lidocaine


90-97846-03

Contents

en	3
bg	8
cs	13
de	17
el	22
es	26
et	30
fi	34
fr	38
hr	42
hu	46
it	50
lt	54
lv	58
nl	62
no/da	67
pl	72
pt	76
ro	80
ru	84
sk	89
sr	93
sv	97
tr	101
Symbols on the packaging	106

Assembly instruction

For further information, see section X

INSTRUCTIONS FOR USE

en

I. COMPOSITION

Cross-linked hyaluronic acid	20 mg/mL
Lidocaine hydrochloride	3 mg/mL
Phosphate buffered saline pH 7	qs ad 1 mL

II. DESCRIPTION

Restylane® Refyne™ is a sterile, biodegradable, transparent gel of non-animal cross-linked hyaluronic acid with the addition of lidocaine hydrochloride 3 mg/mL. The gel is supplied in a prefilled plastic syringe. The contents of the syringe are sterilized using moist heat. The syringe is packaged individually in a blister, with two 30G x 1/2" Ultra thin wall needles. The needles have been sterilized using irradiation. The product is for single use only. To ensure traceability the package includes patient record labels that should be attached to patient records.

III. INTENDED USE

The product is intended to augment the volume of facial tissues. It is recommended to be used for correction of moderate wrinkles, or to redefine the shape of the lips, cheeks or tear troughs. Depending on the area to be treated and the tissue support, the product should be injected in the mid to deep dermis, submucosal layer of the lips, subcutis or supraperiosteal zone. Lidocaine is added to the formulation to diminish the pain resulting from the injection during the treatment.

The product is only intended to be used by authorized personnel in accordance with local legislation, trained in the appropriate injection techniques. Before the first treatment session, it is recommended to contact your local Galderma representative or Restylane distributor for information about training opportunities.

IV. MODE OF ACTION

The product adds volume to the tissue thereby restoring the skin contours of the face to the desired level. The volume and lifting capacity originate from the ability of cross-linked hyaluronic acid to bind water.

V. CONTRAINDICATIONS

- Do not use in patients with a history of hypersensitivity to streptococcal proteins, as the product may contain trace amounts of such material.
- Do not use in patients with severe allergies manifested by a history of anaphylaxis or history of presence of multiple severe allergies.
- Do not use in patients with known hypersensitivity to lidocaine or to amide-type local anaesthetics.

VI. WARNINGS

- Use at specific sites where there is active disease, such as inflammation (skin eruption such as cysts, pimples, rashes or hives), infection or tumours, in or near the intended treatment site should be avoided until the underlying process has been controlled.
- This product must not be injected intramuscularly or intravascularly. Localized superficial necrosis and scarring may occur after injection in or near vessels, such as the nose and glabellar area. It is thought to result from the injury, obstruction, or compromise of blood vessels. Special caution should be taken if the patient has undergone a prior surgical procedure in the planned treatment area. Areas with limited collateral blood flow has an increased risk of ischemia. Aspiration prior to injection is recommended.

- Unintentional introduction of soft tissue fillers into the vasculature in the face may lead to embolisation, occlusion of the vessels, ischemia, necrosis or infarction at the implant site or in the area supplied by the blood vessels affected. Rare but serious adverse events include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral haemorrhage leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if any of the following symptoms occurs, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
- Patients with bleeding disorders or patients using substances that affect platelet function, thrombolytics or anticoagulants may, as with any injection, experience increased bruising or bleeding at the injection site.
- This product should not be mixed with other products prior to injection.

VII. PRECAUTIONS

- This product should only be used by health care practitioners who have appropriate training, experience, and knowledge about the anatomy at and around the site of injection in order to minimize the risks of potential complications (perforation or compression of vessels, nerves and other vulnerable structures).
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Patients with unattainable expectations are not suitable candidates for treatment.
- The safety for use during pregnancy, in breastfeeding females or in patients under 18 years has not been established.
- Injection procedures are associated with a risk of infection. Aseptic technique and standard practice to prevent cross-infections are to be followed.
- Do not use the product if package is opened or damaged, or if the expiry date or lot number is illegible.
- Do not use the product if the contents of the syringe are cloudy.
- This product should be used with caution in patients on immunosuppressive therapy.
- This product should be used with caution in patients with a tendency to form hypertrophic scars or any other healing disorders.
- At each treatment session a maximum dosage of 2 mL per treatment site is recommended.
- Avoid injecting into areas in close proximity to permanent implants, as this could aggravate latent adverse events or interfere with the aesthetic outcome of the treatment. Limited data is available on injecting into an area where a non-permanent implant other than hyaluronic acid has been placed.
- Injection too superficially, or in facial areas with limited soft tissue support or soft tissue cover, or thin skin, such as the periorbital area, may result in contour irregularities and palpable lumps and/or bluish discolouration.
- Injection in the lower periorbital region in patients with pre-existing pigmented dark lower eyelid circles and pre-existing tendency toward oedema formation may be associated with prominent discolouration and excessive swelling due to fluid build-up.

- Patients should avoid excessive sun, UV lamp exposure and extreme temperatures at least until any initial swelling and redness has resolved.
- If laser treatment, chemical peeling or any other procedure based on active dermal response is performed after treatment with this product there is a theoretical risk of eliciting an inflammatory reaction at the implant site. This also applies if the product is administered before the skin has healed completely after such a procedure.
- Injection procedures can lead to reactivation of latent or subclinical herpes viral infections.
- Post inflammatory pigmentation changes may occur after dermal filler injections in people with dark skin (Fitzpatrick Type IV-VI).
- Individual variation and treatment area may affect the bio-degradation of this product, in rare cases product remnants have been detected in tissue when the clinical effect has returned to baseline.
- This product is packaged for single use. Do not resterilize.
- Considerations should be given to the total dose of lidocaine administered if dental block or topical administration of lidocaine is used concurrently. High doses of lidocaine (more than 400 mg) can cause acute toxic reactions manifesting as symptoms affecting the central nervous system and cardiac conduction.
- Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type local anaesthetics e.g., certain anti-arrhythmics, since the systemic toxic effects can be additive.
- Lidocaine should be used with caution in patients with epilepsy, impaired cardiac conduction, severely impaired hepatic function or severe renal dysfunction.

VIII. ADVERSE EVENTS

Anticipated injection-related reactions

Injection-related reactions (including bruising, erythema, itching, swelling, pain and tenderness at the implant site) might occur after treatment. Typically, resolution is spontaneous within one week after injection.

Post marketing adverse event reporting

The following post marketing adverse events have been reported from worldwide sources after treatment with Restylane Refyne (non-exhaustive list). The frequency of reporting is based on the number of estimated treatments performed with the product.

1/1 000 – 1/10 000 Swelling/oedema with immediate onset and onset up to several weeks after treatment.

1/10 000 – 1/100 000 Mass/induration, Erythema, Papules/nodules, Bruising/bleeding, Device ineffective, Discolouration/hyperpigmentation, Pain/tenderness, Infection/abscess including pustule and purulent discharge, Inflammation, Ischemia/necrosis including pallor and vascular occlusion, Eye disorders including dry eye, eye pain and visual impairment such as blurred vision, Other injection site reactions and skin reactions including burning sensation, exfoliation, irritation and warmth, Hypersensitivity/angioedema, Pruritus, Device dislocation.

<1/100 000 Blisters/vesicles, Neurological symptoms including hypoesthesia and paraesthesia, Deformity/asymmetry, Scar/skin atrophy, Rash, Granuloma/foreign body reaction, Acne, Discharge/extravasation, Urticaria, Capillary disorders such as telangiectasia, Dermatitis, Reactivation of herpes infection, Muscle disorders such as muscle twitching, Non-dermatological events including chills, cutaneous contour deformity, dizziness, dyspnoea, headache and sinusitis, Other dermatological events such as chapped lips.

Vascular compromise may occur due to an inadvertent intravascular injection or as a result of vascular compression associated with implantation of any injectable product. This may manifest as blanching, discolouration, necrosis or ulceration at the implant site or in the area supplied by the blood vessels affected; or rarely as ischemic events in other organs due to embolisation. Rare but serious cases of ischemic events associated with temporary or permanent vision impairment, blindness, cerebral ischemia or stroke have been reported following facial aesthetic treatments.

Symptoms of inflammation at the implant site commencing either shortly after injection or after a delay of up to several weeks have been reported. In case of unexplained inflammatory reactions infections should be excluded and treated if necessary since inadequately treated infections may progress into complications such as abscess formation. Treatment using only oral corticosteroids without concurrent antibiotic treatment is not recommended.

The prolonged use of any medication, e.g., corticosteroids or antibiotics in treatment of adverse events has to be carefully assessed, since this may carry a risk for the patient. In case of persistent or recurrent inflammatory symptoms, consider removal of the product by aspiration/drainage, extrusion or enzymatic degradation (use of hyaluronidase has been described in scientific publications). Before any removal procedure is performed, the swelling may be reduced by using e.g. NSAID for 2-7 days or a short course of corticosteroids for less than 7 days, in order to more easily palpate any remaining product.

For patients who have experienced clinically significant reactions, a decision for retreatment should take into consideration the cause and significance of previous reactions.

For reporting of adverse events contact your local Galderma representative or distributor for this product.

IX. PERFORMANCE

In a randomized controlled multicentre study with Restylane Refyne for correction of moderate nasolabial folds, a majority of subjects maintained a clinically significant improvement up to at least 36 weeks. In another randomized controlled multicentre study with Restylane Refyne for correction of moderate to severe nasolabial folds, effectiveness was maintained up to at least 48 weeks after initial injection. In an open-labelled multicentre study to assess performance and tolerance of Restylane Refyne in facial rejuvenation, most subjects were reported with sustained improvement up to at least 6 months in areas including nasolabial folds, marionette lines, and tear troughs. In another open-labelled multicentre study with Restylane Refyne, subjects were assessed for perception of naturalness during contraction following correction of wrinkles and folds in the lower face. Results showed that 98% of subjects had maintained naturalness for at least up to 6 months.

X. ASSEMBLY OF NEEDLE TO SYRINGE

It is important to use a sterile, appropriate needle or blunt cannula. Suitable needles (30G x 1/2" Ultra thin wall) are supplied with the syringe in the blister pack. As an alternative, a blunt thin walled cannula with a recommended size of 27G can be used.

Use surgical gloves, remove the cap from the needle and unscrew the tip cap from the syringe. Hold firmly around the syringe barrel and grasp the needle shield with the other hand. Screw the needle tight onto the syringe by simultaneously pushing and rotating firmly until the needle is completely locked. To ensure proper assembly, minimize the gap between the needle shield and the syringe. See Figure 1.

Remove the needle shield just before injection by pulling it straight out. Do not rotate.

Note! Improper assembly may cause leakage or needle disconnection.

Assembly instruction applies to both needle and blunt cannula.

XI. TREATMENT PROCEDURE

- Inform the patient about the treatment procedure and expected results. Ensure the patient receives appropriate post-treatment care instructions and is made aware of signs and symptoms of potential complications.
- Cleanse the area to be treated with an antiseptic and allow it to dry before injection.
- To avoid breakage of the needle/cannula, do not attempt to bend it before or during treatment. If the needle gets bent, discard it and complete the procedure with a replacement needle.
- Before injecting the product, depress the plunger rod carefully until a small droplet is visible at the tip of the needle/cannula.
- Align the bevel of the needle by turning the syringe on its axis.
- If a blunt cannula is used, an entry point is made in the skin, for example with a sharp needle of appropriate size. During injection, keep the side hole of the cannula facing downwards, away from the skin surface.
- Aspiration is recommended prior to injection in order to reduce the risk of inadvertent injection into a blood vessel.
- Inject the gel slowly by gently pressing down on the plunger rod with the thumb or palm of the hand. Do not apply excessive pressure to the syringe at any time. Presence of scar tissue may impede advancement of the needle/cannula. If resistance is encountered the needle/cannula should be partially withdrawn and repositioned or fully withdrawn and checked for function.
- Choose from a variety of injection techniques, i.e. serial puncture, linear threading or cross-hatching. It is recommended to change needle/cannula for each new treatment site.
- Defects should be fully corrected, but not overcorrected, at each treatment session. Gently massage the treated area after injection.
- If "blanching" of the skin is observed, the injection should be stopped immediately and the whitened area should be massaged gently until it returns to a normal colour before continuing with the injection.
- If the treated area is swollen directly after the injection, an ice pack with adequate protective cloth can be applied on the site for a short period. Ice should be used with caution if the area is still numb from anaesthetic to avoid thermal injury.
- Additional treatments may be necessary to achieve and maintain the desired level of correction.

The syringe, disposable needle/cannula and any unused material must be discarded immediately after the treatment session and must not be reused due to risk for contamination of the unused material and the associated risks including infections. Disposal should be in accordance with accepted medical practice and applicable national, local or institutional guidelines.

XII. SHELF LIFE AND STORAGE

Do not use after the expiry date indicated on package. Store up to 25°C. Protect from freezing and sunlight.

XIII. MANUFACTURER

Q-Med AB, Seminariegatan 21, SE-752 28 Uppsala, Sweden

Phone +46(0)18 474 90 00, Fax +46(0)18 474 90 01

www.galderma.com, e-mail: info.q-med@galderma.com

Restylane and Galderma are registered trademarks.

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Q-Med AB, Seminariegatan 21, SE-752 28 Uppsala, İsveç